

## Microbial Source Tracking Assays Using Methanogen Indicators

**Step 1:** Collect water samples, concentrate bacteria and extra DNA. See Sample Preparation below.

**Step 2:** Determine whether fecal indicator of interest is present by PCR. See PCR Protocols below.

Sewage-specific indicator:

Target organism: *Methanobrevibacter smithii*

Gene target: *nifH*, encodes nitrogenase reductase

Primer name: Mnif

Reference: Ufnar, J.A., S.Y. Wang, J. Christiansen, H. Yampara-Iquise, C. Carson and R.D. Ellender. 2006. Detection of the *nifH* gene of *Methanobrevibacter smithii*: a potential tool to identify sewage pollution in recreational waters. *J. Applied Microbiol.* 101:44-52.

Swine-specific fecal indicator:

Target organism: methanogens

Gene target: *mcrA*, encodes  $\alpha$  subunit of methyl coenzyme M reductase

Primer name: P23-2

Reference: Ufnar, J.A., D.F. Ufnar, S.Y. Wang and R.D. Ellender. 2007. Development of a swine-specific fecal pollution marker based on host differences in methanogen *mcrA* genes. *App. Environ. Microbiol.* 73:5209-5217.

Domesticated ruminant-specific fecal indicator:

Target organism: *Methanobrevibacter ruminantium*

Gene target: *nifH*, encodes nitrogenase reductase

Primer name: Mrnif

Reference: Ufnar, J.A., S.Y. Wang, D.F. Ufnar and R.D. Ellender. 2007. *Methanobrevibacter ruminantium* as an indicator of domesticated ruminant fecal pollution in surface waters. *App. Environ. Microbiol.* (Accepted for publication 8-30-2007).

### Sample Preparation:

1. Collect 500 mL water, place on ice and transport to lab.
2. Remove particulate matter with 3  $\mu$ m filter. Filter 500 mL water sample through a 3  $\mu$ m (47 mm) cellulose acetate filter (Pall Gelman #66387; VWR Cat #28149-634). Collect filtrate in sterile flask
3. Concentrate bacteria with 0.2  $\mu$ m filter. Filter filtrate from previous step using a 0.2  $\mu$ m (47 mm) Supor-200 filter (Pall Gelman #66234; VWR catalog # 28147-979).
4. Place filter in a 100 mL sterile beaker with 5-8 mL sterile PBS and stir using a sterile magnetic stir bar on top of the filter for 5-10 min.
5. Collect resuspended cells in centrifuge tube and pellet cells by centrifugation.
6. Discard supernatant and resuspend cells in 65  $\mu$ L sterile water. The sample may be stored at -20°C or processed immediately for DNA extraction.
7. Extract DNA using MO BIO PowerSoil Kit and quantify DNA concentration.
8. Perform PCR using 10 and 20 ng of DNA.

## PCR Protocols:

### Sewage-specific indicator:

Target organism: *Methanobrevibacter smithii*

Gene target: *nifH*, encodes nitrogenase reductase

Primer name: Mnif

Reference: Ufnar, J.A., S.Y. Wang, J. Christiansen, H. Yampara-Iquise, C. Carson and R.D.

Ellender. 2006. Detection of the *nifH* gene of *Methanobrevibacter smithii*: a potential tool to identify sewage pollution in recreational waters. *J. Applied Microbiol.* 101:44-52.

### PCR:

	Stock Conc	Volume	Final Conc
Water		10.8 uL	
10X Buffer	10X	2 uL	1X
BSA	1 %	2 uL	0.1%
dNTPs	2 mM each	2 uL	0.2 mM each
Template DNA*		1 uL	
Primers	5 uM	2 uL	0.5 uM
Thermal stable DNA pol	5 U/uL	0.2 uL	1 U
Total		20 uL	

\*Template DNA: 10- 20 ng per reaction

Primers: Mnif-342f 5'-AACAGAAAACCCAGTGAAGAG-3'

Mnif-363r 5'-ACGTAAAGGCACTGAAAAACC-3'

### Amplification protocol:

Denaturation: 92°C for 2 min

Amplification, 30 cycles:

92°C for 1 min

55.1°C for 30 sec

72°C for 1 min

Final extension:72°C for 6 min

## PCR Protocols:

### Swine-specific fecal indicator:

Target organism: methanogens

Gene target: *mcrA*, encodes  $\alpha$  subunit of methyl coenzyme M reductase

Primer name: P23-2

Reference: Ufnar, J.A., D.F. Ufnar, S.Y. Wang and R.D. Ellender. 2007. Development of a swine-specific fecal pollution marker based on host differences in methanogen *mcrA* genes. *App. Environ. Microbiol.* 73:5209-5217.

### PCR:

	Stock Conc	Volume	Final Conc
Water		10.8 uL	
10X Buffer	10X	2 uL	1X
BSA	1 %	2 uL	0.1%
dNTPs	2 mM each	2 uL	0.2 mM each
Template DNA*		1 uL	
Primers	5 uM	2 uL	0.5 uM
Thermal stable DNA pol	5 U/uL	0.2 uL	1 U
Total		20 uL	

\*Template DNA: 10 -50 ng per reaction

### Primers:

P23-2f-5'-TCTGCGACACCGGTAGCCATTGA-3'

P23-2r-5'-ATACACTGGCGACATTCTTGAGGATTAC-3'

### Amplification protocol:

Denaturation: 92°C for 2 min

Amplification, 30 cycles:

92°C for 30 sec

60°C for 15 sec

72°C for 30 sec

Final extension: 72°C for 6 min

An Internal Amplification Control (IAC) can be included to guard against false negatives. Use  $10^{-9}$   $\mu$ M in the reaction. The sequence of the IAC is 5'- GAGACCGCAG TCTGCGACAC CGGTAGCCAT TGAACCTGCG ATACCGGTTG CTGCTGCTGC AACTGTTGCC TTGTCCTTGA TGTAGTCGAT TGCATAGTAT GTGTAATCCT CAAGAATGTC GCCAGTGTAT GCTGCTGTAG -3' It can be custom synthesized by oligo vendors such as IDT, Inc.

## PCR Protocols:

### Domesticated ruminant-specific fecal indicator:

Target organism: *Methanobrevibacter ruminantium*

Gene target: *nifH*, encodes nitrogenase reductase

Primer name: Mrnif

Reference: Ufnar, J.A., S.Y. Wang, D.F. Ufnar and R.D. Ellender. 2007. *Methanobrevibacter ruminantium* as an indicator of domesticated ruminant fecal pollution in surface waters. App. Environ. Microbiol. (Accepted for publication 8-30-2007).

### PCR:

	Stock Conc	Volume	Final Conc
Water		10.8 uL	
10X Buffer	10X	2 uL	1X
BSA	1 %	2 uL	0.1%
dNTPs	2 mM each	2 uL	0.2 mM each
Template DNA*		1 uL	
Primers	5 uM	2 uL	0.5 uM
Thermal stable DNA pol	5 U/uL	0.2 uL	1 U
Total		20 uL	

\*Template DNA: 10 - 20 ng per reaction

### Primers:

Mrnif-f-5'-AATATTGCAGCAGCTTACAGTGAA-3'

Mrnif-r-5'-TGAAAATCCTCCGCAGACC-3'

### Amplification protocol:

Denaturation: 92°C for 2 min

Amplification, 30 cycles:

92°C for 30 sec

60°C for 15 sec

72°C for 30 sec

Final extension: 72°C for 6 min

An Internal Amplification Control (IAC) can be included to guard against false negatives. Use 10<sup>-9</sup> μM in the reaction. The sequence of the IAC is 5'- TACAGTAGCT AATATTGCAG CAGCTTACAG TGAAGACAAT AAGAAAGTCA TGGTTATTGG CTGCCTTGAA GAGGAATTTG ATGTAATCTT ATATGATGTT CTTGGAGATG TGGTCTGCGG AGGATTTTCA GTTCCTCTAA – 3'. It can be custom synthesized by oligo vendors such as IDT, Inc.